



MODEL A9070

Shown with Optional 69 Gas Spring Backrest

**Model A9070
Dual Function Split Leg Lift Table**



- Extremely versatile and functional split leg design.
 - Each padded leg rest can be adjusted with knee fully extended and locked into desired angle up to 45° with spring-loaded lock-in lever.
 - Each padded leg rest can be adjusted with hip and knee joints locked into 45° angle with easy to operate hinge mechanism.
- Open shelf storage, accessible from both sides of the table
- Plywood top for added strength.
- 30" x 78" upholstered top with 2" high-density foam.
- Natural Oak, Wild Cherry, or Folkstone Gray laminate.
- Choice of 9 Pro-Form™ vinyl upholstery colors.
- Full subtop for increased strength.
- Available with optional 69 Gas Spring Backrest. (400 lbs. weight capacity)



MODEL	LENGTH	WIDTH	HEIGHT
A9070	78"	30"	31"