

AutoSound™ 7.6 Combo – Automatic Hands-Free Ultrasound & Stimulation Combo

Revolutionary new, patented design allows you to apply ultrasound with a stationary applicator that automatically moves the ultrasound output for a better, more efficient ultrasound treatment. Now combine it with two channels of stimulation for Hands-Free Combo or concurrent ultrasound and stimulation treatments.

- Patented AutoSound design actually moves the ultrasound output for you at the ideal rate of 4cm/second
- Never stops moving the ultrasound output or moves it too fast
- Focus' ultrasound treatment directly over the area – no problems with trying to treat too large an area
- The unique AutoGel Pads, solid conductive gel, that fit perfectly in the Gel Well and make treatment much easier, cleaner, and quicker
- You have three different size straps that come standard with the device to allow the easiest strapping of the applicator to the patient.
- The AutoSound is extremely safe with monitoring circuits plus, a Patient Control Switch for added safety
- Optional manual transducer that provides both a 2cm² and 5cm² soundhead available.

AutoSound Ultrasound Specs

- The AutoSound Hands-Free Applicator is capable of 1MHz, 3MHz, or 1&3MHz Chain that alternates frequency to treat the entire tissue depth.
- Capable of 2.0 W/cm² with 1MHz and 1.5 W/cm² with 3MHz or the 1&3MHz Chain
- Continuous and Pulsed Outputs (100%, 50%, 20%, and 10%)
- AutoSound ERA: 3.5 cm², AutoSound treatment area: 14 cm²
- Timer up to 20 minutes
- 15 programmable ultrasound presets
- Optional Therapy Hammer manual transducer (2cm² and 5cm² on single transducer)
- BNR – Maximum of 5.5 : 1

Stimulation

- Stimulation timer up to 99 minutes
- Two channels of stimulation
- Five simple, preset stimulation treatments that cover 90% of stimulation applications:
 - P1 Quadpolar IFC, Vector, 80-150Hz pulse rate scan
 - P2 Premod IFC, Vector, 2-200Hz pulse rate scan
 - P3 High Volt, 100 Hz fixed pulse rate
 - P4 Premod IFC, 10/10 sec Alternating, 50Hz fixed pulse rate
 - P5 Russian, 10 seconds on / 50 seconds off Surge, 50Hz fixed pulse rate

The AutoSound Combo is completely upgradeable and will periodically have FREE stimulation upgrades to expand the stimulation capabilities!

Combination

- The AutoSound Combo can provide Hands-Free Combination by using the AutoSound applicator to provide ultrasound and as a stimulation electrode.
- Stimulation presets P1, P2, or P3 are available for Combination therapy

Concurrent Stimulation and Ultrasound Treatments

- The AutoSound Combo allows you to provide both a Hands-Free Ultrasound and a separate stimulation treatment at the same time.

Dimensions

12"W x 12"D x 7"H

Weight

Approx 14 lbs