



HCI FITNESS

MSRP: \$599
Call Now!
800-398-9121
360-321-7090

CONTACT US

Hydraulic UBE Table
Adjustable Height UBE Table



Ask About Our UBE & Table Packages!



[Click to Print Hydraulic UBE Table Specifications](#)

Hydraulic UBE Table	
Technical Details	
Hydraulic Adjustment	✓
Adjustment Range	27" - 37"
T-Shaped Table	✓
Table Width	17" / 8.5"

HCI Fitness now offers a custom-made **Hydraulic UBE Table** that is specifically designed and approved for use with the Monark 881e, Monark 871e and HCI PhysioTrainer.

In addition, our **Hydraulic UBE Table** can be used with some other table top ergometers. Refer to the dimensions of your unit, some modification may be needed.

Ask about our special **UBE and Table Packages!**

Product Videos

[PhysioStep LXT - Linear ...](#)

MORE VIDEOS

Table Length	27"
Outside Base Width	30"
Inside Base Width	25.5"
Commercial Casters	✓
Locking Stabilizers	✓
Foot Plate Traction	✓
Mount System	4 Locking Brackets
Product Type	Commercial Grade
Model #	Htable
Warranty	
Frame	Lifetime
Shipping & Assembly	
Shipping	UPS
Assembly Required	Center Support, Hydraulic Cylinder, Table Top
Unit Dimensions	
Compact Dimensions	28" L X 30" W X 27"-37" H
Unit Weight	66 lbs.

Related Products



PhysioTrainer

[ALL PRODUCTS](#)

Testimonials



Customer Reviews

Read product reviews and customer testimonials for HCI Fitness products!

[READ REVIEWS](#)

[PRODUCT BROCHURE](#)

[USER MANUAL](#)

[CONTACTSOURCES](#)