



HCI FITNESS

MSRP: \$4,999

Call Now!

800-398-9121

360-321-7090

CONTACT US

PhysioMax

Independent Bi-Directional Upper Body Ergometer and Recumbent Bike
 all-in-one Total Body Trainer



The **PhysioMax** Independent Bi-Directional Upper Body Ergometer and Recumbent Bike all-in-one Total Body Trainer is the newest innovation in affordable Total Body rehabilitation products with patent pending "independent action."

The **PhysioMax's** independent arm and leg motion allows users to work their upper body (Upper Body Ergometer) & lower body (recumbent bike) independently from each other.

The resistance on the **PhysioMax** is bi-directional & can be adjusted providing versatility for any application. It features a walk-thru design with a removable seat for wheelchair access. The oversize seat swivels to accommodate users with low mobility.

The **PhysioMax** Independent Bi-Directional Upper Body Ergometer and Recumbent Bike all-in-one Total Body Trainer has adjustable hand cranks, pedals, as well as an adjustable UBE arm can be raised or lowered for variable arm positions. Resistance varies from the low "easy start" for deconditioned users and goes up to level 16. The large 9" LCD console has preset programs & displays speed, time, distance, pulse, RPM, Watts,



[Click to Print PhysioMax Specifications](#)

PhysioMax

Technical Details

Swivel Seat	18" Oversized Seat Rotates 360° for Easy Access
Seat Position	Fully Adjustable Forward and Backward 14 Positions
Removable Seat	Seat Removes for Wheelchair Access (ADA)
Resistance System	Electromagnetic

Resistance Levels	1-16
Programs	18 Set Programs (See Computer Details)
Handles	Adjustable Comfort Grip Handles
Pedals	Othopedic Rehab Pedals with Heel Support
Stride	Independent Cycle Motion
Max User Weight	500 lbs.
Power	AC 110V-240V 50/60Hz
Product Type	Commercial Grade
Model #	TBT-1000
Computer Details	
Display	High Definition 9" LCD Display
Backlight	Easy to Read High Contrast Blue
USB Summary	Export Workout Data Including Peaks and Averages
Feedback	Time, Distance, Pulse, Steps Per Minute (SPM), Calories, WATTS, METs, & Resistance Level "1-16"
Programs	18 Programs: Manual, Warm Up, Interval, Ramps, Mountain, Rolling, Climb, Random, Hill, Fat Burn, Cardio, Strength, HRC 65%, HRC 85%, Constant SPM, Constant Watts, Goal Calories, Goal Distance
Heart Rate	Built in Receiver Compatible with Polar and other Manufactures *Chest Belt Not Included
Warranty	

METS & calories.

Product Videos

PhysioMax - Total Body T...

PhysioMax - Total Body T...

MORE VIDEOS

Related Products



PhysioStep RXT

ALL PRODUCTS

Testimonials



Customer Reviews

Read product reviews and customer testimonials for HCI Fitness products!

Frame	Lifetime
Mechanical Parts	3 Years
Electronic Parts	3 Years
Labor	1 Year
Shipping & Assembly	
Shipping	Freight
Assembly Required	Seat, Arms, Display
Unit Dimensions	
Compact Dimensions	60" L X 25" W X 59" H
Unit Weight	350 lbs.

PRODUCT BROCHURE

USER MANUAL

CONTACTSOURCES