

REX RECOVERY EXERCISE X-TRAINER

Combo Pneumatic Compression System

The REX is a four part sequential leg compression system that allows for faster, clinically proven recovery. This fast compression boosts circulation and lymphatic drainage, as well as relieves muscle fatigue. The REX includes a pump and lower extremity garments with air cells. During treatment the pump inflates the air cells to compress muscles, promoting venous blood return and helping prevent deep venous thrombosis. This system is of high quality and durability, ensuring a long lasting experience and uninterrupted performance.

REX

recovery exercise X-trainer



2 treatment modes
3 time settings



Adjustable
pressure
settings



4 Part
Leg Sequential
Compression

FEATURES:

- Two Treatment Modes
- Three Time Settings - 10, 20 and 30 Minutes
- Adjustable Pressure Settings
- Four Part Sequential Compression
- Universal and Large Sizing
- Includes Richmar Backpack to Take Your REX on the Go

ORDER DETAILS:

DVTREX-L	REX DVT System Large
DVTREX-U	REX DVT System Universal
BAG-RCHMR	Richmar Backpack

