

Material Safety Data Sheet



U.S. Department of Labor

May be used to comply with
 OSHA's Hazard Communication Standard,
 29 CFR 1910.1200. Standard must be
 consulted for specific requirements.

Occupational Safety and Health Administration
 (Non-Mandatory Form)
 Form Approved
 OMB No. 1218-0072

IDENTITY <i>(As Used on Label and List)</i> AUTOGEL PAD	Note: Blank spaces are not permitted. If any item is not applicable, or no information is available, the space must be marked to indicate that.
--	---

Section I

Manufacturer's Name RICH-MAR CORPORATION	Emergency Telephone Number 918-543-2222
Address <i>(Number, Street, City, State, and ZIP Code)</i> 15499 E. 590 TH RD	Telephone Number for Information 918.543.2222
INOLA, OK 74036	Date Prepared FEBRUARY 2003
USA	Signature of Preparer <i>(optional)</i>

Section II - Hazard Ingredients/Identity Information

Hazardous Components (Specific Chemical Identity; Common Name(s))	OSHA PEL	ACGIH TLV	Other Limits Recommended	% <i>(optional)</i>
NONE				

Section III - Physical/Chemical Characteristics

Boiling Point	N/A	Specific Gravity (H ₂ O = 1)	1.02
Vapor Pressure (mm Hg.)	N/A	Melting Point	N/D
Vapor Density (AIR = 1)	N/A	Evaporation Rate (Butyl Acetate = 1)	N/A
Solubility in Water SOLUBLE			
Appearance and Odor BUFF COLORED, TRANSLUCENT GEL PAD; LITTLE OR NO ODOR			

Section IV - Fire and Explosion Hazard Data

Flash Point (Method Used) NON-FLAMMABLE	Flammable Limits N/A	LEL N/A	UEL N/A
Extinguishing Media USE APPROPRIATE EXTINGUISHING MEDIA FOR SURROUNDING FIRE			
Special Fire Fighting Procedures NONE			
Unusual Fire and Explosion Hazards NONE			

(Reproduce locally)

OSHA 174, Sept. 1985

Section V - Reactivity Data

Stability	Unstable	Conditions to Avoid	NONE
Incompatibility (<i>Materials to Avoid</i>) NONE			
Hazardous Decomposition or Byproducts N/A			
Hazardous Polymerization	May Occur	Conditions to Avoid	NONE
	Will Not Occur	X	

Section VI - Health Hazard Data

Route(s) of Entry:	Inhalation? NO	Skin? NO	Ingestion? YES
Health Hazards (<i>Acute and Chronic</i>) NONE KNOWN			
Carcinogenicity:	NTP? NO	IARC Monographs? NO	OSHA Regulated? NO
Signs and Symptoms of Exposure NOT APPLICABLE			
Medical Conditions Generally Aggravated by Exposure NONE KNOWN			
Emergency and First Aid Procedures N/A			

Section VII - Precautions for Safe Handling and Use

Steps to Be Taken in Case Material is Released or Spilled SLIPPERY IF SPILLED ON FLOOR. MOP AREA WITH WATER
Waste Disposal Method FOLLOW ALL FEDERAL, STATE, AND LOCAL REGULATIONS FOR NON-HAZARDOUS WASTE DISPOSAL.
Precautions to Be taken in Handling and Storing DO NOT STORE BELOW 45° f (5°C) OR ABOVE 135°F (57°C)
Other Precautions NONE

Section VIII - Control Measures

Respiratory Protection (<i>Specify Type</i>) NOT APPLICABLE			
Ventilation	Local Exhaust N/A	Special	N/A
	Mechanical (<i>General</i>) N/A	Other	N/A
Protective Gloves N/A	Eye Protection	N/A	
Other Protective Clothing or Equipment N/A			
Work/Hygienic Practices N/A			

ND - NOT DETERMINED

N/A - NOT APPLICABLE OR NOT AVAILABLE

PAGE 2

* U.S.G.P.O.: 1986 - 491 - 529/45775