

Theramini 3C (Clinical) Combo Unit

Clinical Screen Controlled 2 Channel Stimulator and Ultrasound Combination Unit.

- Theramini 3C has 24 nameable, programmable QuikSets for customized treatments that can be entered and even named to meet the facility's needs allowing for quick and easy treatment setup by the whole staff
- **Ultrasound** has 2cm² and 5cm² soundheads on single, 4-in-1 "Therapy Hammer" transducer. Both 2cm² and 5cm² soundheads can output both 1 & 3 MHz frequencies.
- **Stimulator** has two stimulation channels that treat two site with independent intensity levels at the same time.
- Both stimulator and ultrasound can run separately at the same time.
- Upgradeable microprocessor controlled 6 Stimulation Waveforms
- Waveforms available on *both channels* include:
 - **Quadpolar Interferential** (5,000 Hz carrier frequency with Cont & Chain Modes)
 - **Premodulated Bipolar Interferential** (5000 Hz carrier frequency with Cont, Surge, Alt and Chain Modes)
 - **Monophasic** (High Volt Pulsed Current twin square wave with interphase interval adjustment, with Cont, Surge, Alt and Chain Modes)
 - **Biphasic** (Symmetric square wave with phase duration and interphase interval adjustments, with Cont, Surge, Alt and Chain Modes)
 - **Russian** (2500 Hz carrier frequency with Cont, Surge and Alt Modes)
 - **Microcurrent** (Biphasic output, 1 to 1000 Hz pulse rate with Cont and Chain Modes)
- Combination ultrasound and stimulation treatments can be used together or independently from each at the same time.
- Patient safety features include patient isolating safety loads to protect the patient from electrical shock due to power surges or equipment malfunction. Also, a built-in lead cord check
- **3C** - is enclosed in an all metal clinical cabinet
- Comes with 2 each patient lead cords, one package of electrodes and 16 oz bottle of Rich-Mar Aloe Sound Ultrasound Lotion
- One year warranty
- 19" W x 10"D x 2.755" H – 13 lbs. Available in 120 VAC, 60Hz or 220 VAC, 50 Hz

Optional Accessories

Theramini Accessory Kit

- 4 each 4" round, reuseable carbon electrodes
- 4 each sponges for electrodes
- 4 each elastic velcro straps for electrodes

Thermini 3C Carrying Case

Soft, padded nylon case available for added portability. Case has shoulder strap, handles and numerous pockets.

Two Year Premium Extended Warranty

Provides two year extension of warranty that covers all parts and labor plus will provide a Loaner unit via 2nd Day Air if unit breaks. Will also arrange for pick-up and transport of broken unit back to factory for factory repair. Will then ship repaired unit back to customer and arrange for pick up and transport of loaner unit. No cost for labor, parts, loaner or shipping for the entire warranty period. Also, includes FREE upgrades of unit in the warranty period. Premium Extended Warranty is also transferable.