

Therasound 3.5 Ultrasound – Affordable, Advanced & Portable – 1 & 3 MHz Frequencies with patented “Therapy Hammer” Transducer

The Therasound 3.5 offers fantastic flexibility in a lightweight ultrasound that will work equally well in the clinic or on the road.

- Rich-Mar’s patented “Therapy Hammer” transducer that provides single transducer with both a 2cm² and a 5cm² soundhead.
- Both the 2cm² and 5cm² can output both 1 & 3 MHz frequencies.
- The Therasound 3.5 is an affordable, advanced and portable ultrasound that offers both continuous and pulsed outputs.
- Additionally the 3.5 allows you adjust at any time during the treatment:
 - time
 - frequency
 - transducer size
 - intensity
 - output readout.
- There is also a convenient pause feature.
- The clear, bright digital readouts and simple panel make the Therasound 3.5 easy to use.
- The Therasound 3.5 comes with a standard 2 year warranty which includes the cable and transducer (most other manufacturers will only provide one year of warranty on the cable and transducer because that is the part most susceptible to breaking).

Optional Warranty Upgrade

You can even upgrade to a 3 Year Warranty that also includes a 2nd Day Air Loaner FREE for the entire warranty period - A great deal at only \$99 for the warranty upgrade.

Therasound 3.5 Specs

Output: 0-2.0 W/cm² (0-10 Watts)

Power Input: 110 VAC / 60 Hz, (220 VAC / 50 Hz available)

Transducers: Therapy Hammer 2cm² and 5cm² Effective Radiating Areas (ERA)

Frequency: 1 & 3 MHz available on both size soundheads

Beam Non-uniformity Ratio (BNR): 5.5 : 1 maximum

Dimensions: 7.5" W x 7.5" D x 4"H

Weight: 5.6 lbs (2.54 kg)

Duty Cycles: 100%, 50%, 20%, 10%

Time: 0-30 minutes

Warranty: 2 Years on entire unit, including cable and transducer

Fuse: 1 amp

Line Leakage: Less than 50 mA

Therasound 3.5P Dimensions

14"W x 10.5"D x 6.5"H – 12 lbs Hard plastic case that has the ultrasound built into it.